INSTRUCTIONS

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MM6-ILW/ILG/ILPS-2

LED ILLUMINATOR UNIT

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This instruction manual is for the Olympus LED Illuminator unit model MM6-ILW/ILG/ILPS-2. To ensure the safety, obtain optimum performance and to familiarize yourself fully with the use of this attachment, we recommend that you study this manual thoroughly before operating the unit. Retain this instruction manual in an easily accessible place near the work desk for future reference.

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A X 7 9 6 9

This product is applied with the requirements of standard IEC/EN61326-1 concerning electromagnetic compatibility.

- Emission Class A, applied to industrial environment requirements.

- Immunity Applied to industrial environment requirements.

Some interference may occur if this product is used in domestic location.

In accordance with European Directive on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately.

Refer to your local Olympus distributor in EU for return and/or collection systems available in your country.

NOTE: This product has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the product is operated in a commercial environment. This product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this product in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

FCC WARNING: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the product.

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IMPORTANT

The LED Illuminator Unit has been developed to replace the Light Guide Light Source (MM6-LS22) for use with Olympus Measuring Microscopes (STM6, STM6-LM).

The LED Illuminator Unit features a compact size, low power consumption and non-necessity of lamp bulb replacement.

As the instruction manual for the STM6 and STM6-LM refers to the previous light guide light source when it mentions the light source, please use this manual in place of the existing description on the light source.

Safety Precautions

1. Always use the power cord provided by Olympus.
2. Always ground the equipment.
3. The power supply unit incorporates high-voltage parts inside. Do not open the cover of this unit to prevent electric shock.

Notes on the LED (Light Emitting Diode) Safety

It is discommended to stare at the light directly from the LED light source composed in this product for long hours as this could damage your eye. Use the equipment in accordance with the regulations for Workers Health and Safety, if any.

Safety Symbols

The following symbols are found on this unit. Study the meaning of the symbols and always use the equipment in the safest possible manner.

Symbol	Explanation
	Before use, carefully read the instruction manual. Improper handling could result in injury to the user and/or damage to the equipment.
	Indicates that the main switch is ON.
	Indicates that the main switch is OFF.

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Caution

If the unit is not used properly as instructed in this manual, the safety can no longer warranted and the unit may also fail. Be sure to use the unit as instructed in this manual.

The caution messages in this manual carry the following symbol.

▲ Caution for prevention of injury of the user or damage to the equipment (including properties in the surroundings).

Before disposing of this product, be sure to follow the regulations and rules of your local government.

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1 NOMENCLATURE


Power Supply Unit (MM6-ILPS-2)

- ① Main switch
I : ON. O : OFF.
- ② Brightness control
Turn clockwise to increase the light intensity.
△ : Transmitted light illumination
▽ : Reflected light illumination

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
2 ASSEMBLY

▲ Cables and cords are vulnerable to bend or twist. Do not apply excessive force to them.

▲ Always use the power cord provided by Olympus. If no power cord is provided, please select the power cord by referring to the section "PROPER SELECTION OF THE POWER SUPPLY CORD" at the end of this instruction manual. If the proper power cord is not used, Olympus can no longer warrant the electrical safety performance of the equipment.

1 Installing the LED Light Emitter Block (Fig. 1)

1. When installing the White LED Unit (MM6-ILW) ①, insert it into the vertical illuminator unit and clamp with the clamping knob ②.
2. Fit the cable of the White LED Unit in the holder ③, giving it a certain slack so that it does not hinder the vertical movement of the Z-axis.
3. When installing the Green LED Unit (MM6-ILG) ④, insert it into the lower part of the main body and clamp with the clamping knob ⑤.


Fig. 1

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

2

Connecting the Power Supply Unit (Figs. 2 & 3)

▲Ensure that the main switch of the Power Supply Unit is set to “○” (OFF) before proceeding to the connection.

1. Insert the connector ① firmly into the OUTPUT connector ②.
 - White LED Unit: Insert into the connector marked .
 - Green LED Unit: Insert into the connector marked .
2. Connect the power cord ④ into the power inlet receptacle ③, and connect the power plug ⑤ of the power cord into the power outlet ⑥.

▲Connect the provided power cord correctly and ensure that the grounding terminal of the power supply and that of the 3-conductor wall outlet are properly connected. If the equipment is not grounded, Olympus can no longer warrant the electrical safety performance of the equipment.


Fig. 2

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Power Supply Unit

Model	MM6-ILPS-2
System	Constant-current lighting, 2 circuits
Output current variation	0 to 630 mA
Rated supply voltage	100-240 V AC 0.2-0.1 A, 50/60 Hz
Power consumption	10 W
Dimensions	95(W) x 132(D) x 55(H) mm (excluding projections)
Weight	550 grams

Common Specifications to MM6-ILW/ ILG/ILPS-2

Operating environment

- Indoor use.
- Altitude: Max. 2000 meters.
- Operating temperature/humidity: 5° to 35°C (50° to 95°F) 25% to 85%.
- Storage temperature/humidity: -20° to 60°C (-4° to 140°F) 20% to 85%.
- Supply voltage fluctuations: \pm 10%.
- Pollution degree: 2 (in accordance with IEC60664-1).
- Installation (overvoltage) category: II (in accordance with IEC60664-1)

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■ PROPER SELECTION OF THE POWER SUPPLY CORD

If no power supply cord is provided, please select the proper power supply cord for the equipment by referring to "Specifications" and "Certified Cord" below:

CAUTION: In case you use a non-approved power supply cord for Olympus products, Olympus can no longer warrant the electrical safety of the equipment.

Specifications

Voltage Rating	125V AC (for 100-120V AC area) or, 250V AC (for 220-240V AC area)
Current Rating	6A minimum
Temperature Rating	60°C minimum
Length	3.05 m maximum
Fittings Configuration	Grounding type attachment plug cap. Opposite terminates in molded-on IEC configuration appliance coupling.

Table 1 Certified Cord

A power supply cord should be certified by one of the agencies listed in Table 1, or comprised of cordage marked with an agency marking per Table 1 or marked per Table 2. The fittings are to be marked with at least one of agencies listed in Table 1. In case you are unable to buy locally in your country the power supply cord which is approved by one of the agencies mentioned in Table 1, please use replacements approved by any other equivalent and authorized agencies in your country.

Country	Agency	Certification Mark	Country	Agency	Certification Mark
Argentina	IRAM		Italy	IMQ	
Australia	SAA		Japan	JET	
Austria	ÖVE		Netherlands	KEMA	
Belgium	CEBEC		Norway	NEMKO	
Canada	CSA		Spain	AEE	
Denmark	DEMKO		Sweden	SEMKO	
Finland	FEI		Switzerland	SEV	
France	UTE		United Kingdom	ASTA BSI	
Germany	VDE		U.S.A.	UL	
Ireland	NSAI				

Table 2 HAR Flexible Cord

APPROVAL ORGANIZATIONS AND CORDAGE HARMONIZATION MARKING METHODS

Approval Organization	Printed or Embossed Harmonization Marking (May be located on jacket or insulation of internal wiring)		Alternative Marking Utilizing Black-Red-Yellow Thread (Length of color section in mm)		
			Black	Red	Yellow
Comite Electrotechnique Belge (CEBEC)	CEBEC	⌊HAR⌋	10	30	10
Verband Deutscher Elektrotechniker (VDE) e.V. Prüfstelle	⌊VDE⌋	⌊HAR⌋	30	10	10
Union Technique de l'Electricite' (UTE)	USE	⌊HAR⌋	30	10	30
Istituto Italiano del Marchio di Qualita' (IMQ)	IEMMEQU	⌊HAR⌋	10	30	50
British Approvals Service for Electric Cables (BASEC)	BASEC	⌊HAR⌋	10	10	30
N.V. KEMA	KEMA-KEUR	⌊HAR⌋	10	30	30
SEMKO AB Svenska Elektriska Materielkontrollanstalter	SEMKO	⌊HAR⌋	10	10	50
Österreichischer Verband für Elektrotechnik (OVE)	⌊OVE⌋	⌊HAR⌋	30	10	50
Danmarks Elektriske Materialkontroll (DEMKO)	⌊DEMKO⌋	⌊HAR⌋	30	10	30
National Standards Authority of Ireland (NSAI)	⌊NSAI⌋	⌊HAR⌋	30	30	50
Norges Elektriske Materielkontroll (NEMKO)	NEMKO	⌊HAR⌋	10	10	70
Asociacion Electrotecnica Y Electronica Espanola (AEE)	⌊UNED⌋	⌊HAR⌋	30	10	70
Hellenic Organization for Standardization (ELOT)	ELOT	⌊HAR⌋	30	30	70
Instituto Portages da Qualidade (IPQ)	np	⌊HAR⌋	10	10	90
Schweizerischer Elektro Technischer Verein (SEV)	SEV	⌊HAR⌋	10	30	90
Elektriska Inspektoratet	SETI	⌊HAR⌋	10	30	90

Underwriters Laboratories Inc. (UL)
Canadian Standards Association (CSA)

SV, SVT, SJ or SJT, 3 X 18AWG
SV, SVT, SJ or SJT, 3 X 18AWG

Manufactured by

EVIDENT CORPORATION

6666 Inatomi, Tatsuno-machi, Kamiina-gun, Nagano 399-0495, Japan

Distributed by

EVIDENT EUROPE GmbH

Caffamacherreihe 8-10, 20355 Hamburg, Germany

EVIDENT EUROPE GmbH UK Branch

Part 2nd Floor Part A, Endeavour House, Coopers End Road, Stansted CM24 1AL, U.K.

EVIDENT SCIENTIFIC, INC.

48 Woerd Ave Waltham, MA 02453, U.S.A.

EVIDENT AUSTRALIA PTY LTD

97 Waterloo Road, Macquarie Park, NSW 2113, Australia

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EVIDENT Customer Information Center

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