

Get More Out of Your Digital Microscope



# Intelligent Innovation

Fast failure analysis with guaranteed accuracy and repeatability



## Macro to Micro Versatility

- ▶ Large selection of lenses to find the best magnification, resolution, and working distance for your sample
- ▶ Coded free-angle observation system



3 - 8



## Multiple Observations with a Single Click

- ▶ Change lenses and observation method quickly by pushing a button
- ▶ All observation methods are available at all magnifications



9 - 14



## Be Confident in Your Results with Guaranteed Accuracy and Precision


- ▶ Accurate measurements with a telecentric optical system
- ▶ Both accuracy and repeatability are guaranteed at all magnifications



15 - 18



# Macro to Micro Versatility



The microscope's 20X to 7000X magnification range enables you to conduct high-level, low-magnification overview observations and seamlessly zoom down to the micron level for detailed analysis.

The depth of field and a long working distance give you the flexibility to inspect larger samples, while the free-angle observation system enables you to image your sample from many directions.

# The Tool to Solve Your Challenges

## Rough Inspection and Micron-Level Analysis with One System

In the past, both high-magnification and low-magnification microscopes were needed to complete an inspection. Switching your samples between microscopes took time and required many setting adjustments.




- Better objectives deliver better resolution
- Long working distance
- Deep depth of focus
- Quick and easy lens replacement

**DSX1000**

**Complete your inspection with one easy-to-use system.**

## High-Resolution Images at High Magnification

When inspecting uneven samples, it is important to maintain a safe distance between the lens and sample to keep from damaging it. To see details, you need to increase the magnification, but this typically results in worse resolution.



**DSX1000**

**High-quality images at high magnification with advanced optics.**

## Minimize the Chance of Crashing into Your Sample

If the distance between your sample and the lens is too small, the objective can crash into the sample during analysis, potentially damaging it.



**DSX1000**

**Observe uneven samples without bumping into them.**

# Choose the Best Lens for Your Analysis


Our lineup of 17 objective lenses, including super long working distance and high numerical aperture options, provides the flexibility to obtain a wide range of images.



For more information on our lenses, see pages 27 and 28.


## See the Whole Picture: 20X to 7000X Magnification Range

Seamlessly change magnification from high-level analysis to detailed observation by pushing a button.



## Minimize the Chance of Crashing into Your Sample

The DSX1000 system offers a wide depth of field and a long working distance, so you can observe uneven samples with less chance of causing damage.



SXLOB series

## High Resolution and a Long Working Distance in One Objective

Objectives combining high resolution and long working distance enable you to analyze large, uneven samples, such as automobile and machines parts, that were difficult to inspect in the past using an optical microscope.



XLOB series

## Exceptional Resolution with a 0.95 Numerical Aperture

The DSX1000 digital microscope enjoys the full benefits of optical microscope lenses. Their chromatic aberration correction enables you to see the fine details in your sample.




UIS2 series

## See Your Sample from Many Angles



### Oblique observation ( $\pm 90^\circ$ )

The eucentric optical design maintains a good visual field when tilted or when the stage is rotated, enabling you to observe your sample from many angles. This flexibility frees you from only having the option to observe your samples directly from above, helping you spot hard-to-see defects.



### Rotational observation ( $\pm 90^\circ$ )


The stage rotates 90 degrees for even more flexibility in how you view your sample.





## Always Know Your Angles

The system automatically tracks the inclination and angle rotation information for each image.




Movement of rotating stage

Inclination angle sensor



Inclination angle: 45°

Micro map function



You will not get lost during an inspection. The system displays the area you are currently observing within the whole picture, even in zoom mode.

Easy-to-Use Console



XY stage and Z-drive control with joystick

Quickly move the zoom head by adjusting the fine focus knob

# Multiple Observations with a Single Click

## Sliding Nosepiece



## Six Observation Methods



## Console




The DSX1000 microscope offers flexibility to make your inspection workflow faster and easier. Changing observation is as simple as turning a dial while switching between Six different observation methods requires only the push of a button.

## Instant Switching Saves Time


Conventional systems may only offer one or two observation methods, limiting what you can see in your sample. The DSX1000 microscope offers various observation methods from which you can choose the one that works best for your job.

Supported Observation Methods for Conventional Digital Microscopes

	Observation method A	Observation method B	Observation method C
Lens magnification A	Unsupported	Unsupported	Supported
Lens magnification B	Unsupported	Unsupported	Supported
Lens magnification C	Supported	Conditionally supported	Conditionally supported



Generally, replacing lenses on an optical microscope is cumbersome and some illumination methods may not be supported.



**DSX1000**

**Choose from 6 observation methods, and switch between them with a single click.**

## Change Magnification Quickly and Easily

With some digital microscopes, you need to replace the object lens to adjust the magnification. This can be a slow process, potentially requiring you to remove the camera cable each time and also must be restart the software. During this process, you might lose your view on the object forcing you to spend time navigating back to the correct spot.

The DSX1000 enables you to easily and quickly change magnification from the macro to the micro range, minimizing the chance of losing the target object.

### Quick Magnification Change with a Sliding Nosepiece

You can attach two objective lenses to the head at the same time and quickly change the magnification just by sliding the lens.

### Instantly Switch the Lens Attachment

Objective lenses can be quickly switched you to find the best magnification for your inspection. When the lens is replaced, magnification and visual field information will automatically be updated.





### Fast Motorized Optical Zoom

Optically zoom in and out by turning the console dial. The optical zoom head covers wide range of magnifications with a single objective. It's fully motorized, helping you to eliminate common errors that may occur when manually setting the zoom.




Knob dial




## Change Observation and Illumination Methods by Pushing a Button

With some microscopes, illumination methods are depend on your choice of lens, and changing illumination can be time consuming. The DSX1000 system makes this process fast, simple, and easy—just push a button.




## Simple Light Adjustment with Knob Dial

\*Lighting is adjusted differently depending on the observation method.



## Easy-to-Use Console

The multifunction console helps you do your work quickly. For example, you can easily capture 2D or 3D images or move the XYZ stage with one click.





## Integrated Observation Methods

Easily switch between brightfield (BF), oblique, darkfield (DF), MIX (BF and DF), simple polarization (PO), differential interference contrast (DIC), and contrast enhancement observation functions. This flexibility enables you to handle almost any microscope inspection task.

### MIX (BF+DF)

**Light comes from a ring around the lens**

Easily detect scratches and defects that can be hard to find with a conventional microscope by combining the detection capabilities of darkfield (DF) to the visibility of brightfield (BF).



BF


DF

MIX

### BF (Brightfield)

**Good for flat samples**



On a mirrored surface, scratches look dark against the surface, helping them stand out.



### PO (Polarization)

**Designed for polarizing samples**



By orthogonally laying out two polarization filters, this method enables you to see the contrast and color according to your sample's polarization property.



### OBQ (Oblique)

**Enhance your surface's unevenness**


Use this method to enhance a surface's unevenness by shining the light from only one direction. This method is ideal for uneven or corrugated samples and cutting traces.



### DIC (Differential Interference Contrast)

**Visualize unevenness, foreign particles, scratches, and other defects at the nano level**



This method enables you to visualize surface unevenness at the nano level. It's ideal for inspecting wafers, film, LCD ACF, and glass surfaces.



### DF (Darkfield)

**Best for detecting scratches and similar defects**



Scattering or reflected light is obliquely irradiated on the sample's surface, highlighting dust, scratches, and other objects. Dust and scratches appear bright in the visual field.



### Contrast UP


**Emphasize your sample's contours**

This method enhances the contrast by narrowing the optical element's aperture stop, enabling you to see sharp, vivid images. The bright parts look brighter, while the dark parts look darker.



## Minimize Glare

The adaptor diffuses lighting to help eliminate glare and darkened slopes on a samples like a cylindrical metal surface.




Without adaptor

Within adaptor

## Eliminate Reflections

When observing a film's surface or an object through a transparent medium, such as glass, part of the surface can look very bright. An optical polarization plate is used with the adaptor to eliminate glare.




Without adaptor


Within adaptor

## High-Resolution Images at High Speed

The microscope's advanced algorithms enable you to quickly capture 3D images by pushing a button.



Fully focused image 3D image



Obtaining a high visual field image by connecting images

## Panoramic Images with Automatic Stitching

Capture 3D images over a wide area in panorama view. Stitch together a series of in-focus images to see your sample beyond the microscope's field of view.



Panoramic image

# Be Confident in Your Results with Guaranteed\* Accuracy and Precision



The microscope's telecentric optical system enables you to obtain very precise measurements while the guaranteed accuracy and precision bring confident in your results.

\*To guarantee XY accuracy, calibration work must be undertaken by an Olympus service technician




# Guaranteed Measurement Precision

---


## Be Confident in Your Measurements

The precision of many general digital microscopes and optical microscopes is not guaranteed.

Conventional manual measurement



**DSX1000** with measurement accuracy



**DSX1000**

You can be confident in your measurement results with the guaranteed measurement precision.

## On-Site Calibration

Even if the measurement precision of your microscope was guaranteed at the time of factory shipment, those results can be changed once installed.

Conventionally there is no calibration certificate



**DSX1000** with calibration certificate



**DSX1000**

Reliable measurement with on-site calibration.

## High-Precision Measurement

When imaging tall samples with conventional microscope, you may suffer from a convergence effect where the size of the object can look different depending on the point of focus. This effect makes it difficult to take accurate measurements. The DSX1000 system's telecentric optics eliminate this effect to achieve better measurement accuracy. When you need high-precision measurement, the DSX1000 is your choice.

Conventional digital microscope (non-telecentric optical system)	<b>DSX1000</b> (telecentric optical system)	
The size is different between the right and left edges in one visual field.		The size is the same between the right and left edges in one visual field.

### What is a telecentric optical system?

Telecentric lenses have the same brightness at the center and edge of the visual field. Even if the sample moves vertically by adjusting the focus, the image size (magnification) does not change with telecentric lenses. This optical system enables you to capture an image of an entire sample faced up, which increases measurement precision.


Non-telecentric optical system	<b>Telecentric optical system</b>		
When measuring the distance between two points in the images above and below focus, results can differ.	The measurement result is the same between the images above and below focus.		
	<p>Normal lens</p> <p>With a normal lens, the target surface can be partially hidden by unevenness.</p>	<p>Telecentric lens</p> <p>With a telecentric lens, the target surface is not hidden by unevenness.</p>	
	<p>The images are different in size.</p>	<p>The image size is the same.</p>	
Below focus			Below focus

## Guaranteed Accuracy and Repeatability

Measurement accuracy and repeatability are guaranteed at all magnifications, so you can be confident with your measurement results.

Measurement object: 1.00 mm standard scale

Measurement count	Measurement result
1	1.0 mm
2	1.02 mm
3	0.99 mm
4	1.01 mm
5	1.0 mm
6	1.0 mm
7	0.99 mm
Measurement count	Average value
7	1.00 mm




- To issue certificates, calibration work must be undertaken by Olympus's dedicated service staff.
- Olympus issues the calibration certificate authenticated by the ILAC-MRA calibration accreditation agencies.

## Guaranteed Measurement Performance at Your Working Environment


When you purchase DSX1000 system, calibration will be done by a technician at your site to guarantee the same level of precision as it was shipped from the factory.

A variety of certifications



## Keep Your Measurements Precise

To further reduce fluctuation in measurement precision, the objective lenses and zoom ratios needs to be calibrated. Normally, this is a time-consuming process, but calibration setting can be done quickly and easily by the auto calibration feature.




Calibration sample

# Powerful Functionality Delivers Exceptional Value



## Smooth Live Imaging with a Fast 60 fps Frame Rate

With the same technology used in high-end, digital single-lens reflex cameras, the DSX1000 offers smooth imaging at 60 frames-per-second (fps) frame rate. Your images remain sharp even if you move the sample.



## High-Resolution Imaging for High Color Reproduction

You can obtain high-resolution images with exceptional color reproduction and a small file size with the camera's built-in 3CMOS mode.




The DSX1000 system can achieve the same image quality as a three-plate camera by capturing images successively after shifting the sensor's position.

## Sharp Low-Magnification Images without Flare

The microscope's advanced optical technology eliminates lens flares that are common at low magnification, resulting in sharp images.


## Preview Images from 6 Observation Methods

Instantly display sample images captured with 6 different observation methods by a single click. Choose the image that works best for your sample, and the setting will automatically be configured to make the best out of that observation method.




## Retrieve Optimal Observation Condition

When you take an image, it contains information about the conditions when captured. You can recall these conditions by clicking on the image, making it easy to observe with the same conditions and settings.





## Minimize Halation

The HDR function combines multiple images captured at different exposures to show the fine structures in bright and dark areas while removing halation and glare from reflective samples.




## Wide Variety of Measurements

The system not only supports measurement of 2D properties like line width, surface area, angle, and diameter, but can also measure height, volume, cross-sectional area, and other properties required for 3D measurement.



## Surface Roughness Measurement

You can easily see the picture of surface condition by performing line and surface roughness measurement quantitatively, using Ra and Rz parameters.



解析パラメータ			
Sq	21.104 (μm)	Sk	0.531
Skw	1.996	Sp	46.136 (μm)
Sv	28.662 (μm)	Sz	74.798 (μm)
Sa	18.311 (μm)		

## Powerful Image Analysis Software

OLYMPUS Stream image analysis software facilitates specialized analysis, such as granularity measurement. OLS5000-BWS software is also available to make your inspection, from data acquisition to reporting, more efficient.

## Flexible, One-Click Reporting

Instantly report your results in the format you like. This tool supports Excel, PDF, and RTF formats in addition to DSX-specific formats. You can also customize your report according to your desired format.

## Automotive


Observe foreign substances on the surface of painted car bodies to identify the source of the contamination




A painted car body

### Solution

See clear object details at the same magnification.




The sample is blurry when using a conventional objective (1700X)



Detect foreign substances by clearly visualizing even the air bubbles around them (DSX1000, 1700X)


Observe the cross-section of radiator fins to find welding defects




Radiator fin cross-section

### Solution

With some systems, it can be time-confusing to choose the best observation method for your sample. With the DSX1000 system, select your observation method simply by pushing a button.



The sample is blurry when using a conventional objective (1700X)



DSX1000 Polized observation (300X), clearly visualize the peeling of the weld

# Metal


## Observe a fractured metal surface to analyze the damage cause



A fractured metal surface

### Solution

Observation of a wide area can be done under high magnification with image stitching, but with some conventional system, the boundaries of stitched images remain visible. The improved stitching algorithm of DSX1000 provides clear images without visible boundaries.



A 2 × 2 stitched image (1000X)


## Observe Large Objects without Glare



Drill tip

### Solution

Brightness can easily be inequable when observing large three-dimensional objects, which makes it difficult to see the entire view of a sample. Obtain a clear, glare-free overview of a large object with the DSX1000 microscope.




Using a conventional objective (24X), inequable brightness makes it hard to see the damaged area



Easy observation (24X) of damaged area with DSX1000 thanks to the flat lighting

# Electronics


## Measurement of integrated circuit (IC) chipping to determine the cause of failure



IC wafer before dicing

### Solution

Not every digital microscope guarantees the accuracy and repeatability of measurements at all magnifications. You can rely on the DSX1000 measurement results with guaranteed measurement accuracy and repeatability.



Differential interference contrast (DIC) image (2500X), the edge of chipping is clearly visible


## Multilayer capacitor surface inspection of defects and external dimension measurement



IC substrate

### Solution

Reflection between the capacitor and dielectric makes it challenging for a conventional digital microscope to observe the entire surface. Selecting the appropriate observation method can instantly be done with DSX1000 to find the best image.




Brightfield observation (1500X), surface observation and external dimension measurement can be done at the same time



## Other analysis applications

### Analyze the characteristics and defects in the cross-section of metallic materials




A polished sample


#### Solution

The DSX1000 system with OLYMPUS Stream software is capable of acquiring a fully-focused image of the entire sample regardless of unevenness or tilt on the polished surface.

This eliminates the need for re-polishing, leading to a reduction of effort and time.




Using a conventional objective (100X), the sample is only partially in focus



DSX1000 objective (100X), the entire sample is fully focused regardless of irregularities


### Analyze the glass fibers and resin in the cross-section of a printed wiring board's glass epoxy substrate











Printed wiring board

#### Solution


The glass epoxy's substrate is uneven due to etching, making it difficult for the microscope to clearly focus. The DSX1000 objectives' depth of focus and resolution achieve clear images across the entire cross-section.




Darkfield (700X), individual glass fibers can be clearly observed

Model		Entry model	Tilt model	High Resolution model	High End model	
						
						
Feature / Customer benefit		Basic functionality and easy to operate	Preferred for analyzing irregularly shaped samples	High-resolution images for advanced analysis	Analyze a wide variety of sample types using multiple observation methods	
Standard equipment	Microscope motorized zoom head	Universal zoom head *DIC : Differential interference contrast *Depth of focus up *High resolution 3CMOS mode	—		●	●
		Standard zoom head	●	●	—	
		Observation method BF :Brightfield DF :Darkfield OB :Oblique MIX :MIX POL :Polarized light	●	●	●	●
	Microscope frame	Tilting frame (±90°)	—	●	—	●
		Upright frame	●	—	●	—
	Stage	Motorized XY stage with rotation (±90°)	—		●	●
		Motorized XY stage	—	●	●	—
		Manual XY stage	●	—		—
	Console		—	●	●	●
	Objectives*	Super long working distance objective lens	*Please refer to the objectives lens lineup on page 27 – 28			
Long working distance objective lens						
UIS2 objective lens						
Others	Application software	●	●	●	●	
	Calibration sample	●	●	●	●	
	Controller PC/Display monitor	●	●	●	●	
Option	Transmitted lighting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Adaptor	Diffusion adaptor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Eliminate reflection adaptor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Software	Edge detecting measurement	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Particle analysis	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Others	Objectives strage case	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

● : Standard □ : Option



# System diagram



# Objective Lenses

## Super long working distance objective lens

- Provides a long working distance between the lens and sample



## High-resolution, long working distance objective lens

- Delivers both high resolution and a long working distance




## High-performance, high NA objective lens

- Delivers high performance at nano scale



Magnification on monitor 20X 40X 100X 200X

Objective lens model	20X	40X	100X	200X
DSX10-SXLOB1X	20 – 140X			
DSX10-SXLOB3X		42 – 420X		
DSX10-SXLOB10X				
DSX10-XLOB3X		42 – 420X		
DSX10-XLOB10X				
DSX10-XLOB20X				
DSX10-XLOB40X				
MPLFLN1.25X	22.5 – 175X			
MPLFLN2.5X		37.5 – 350X		
MPLFLN5XBDP			70 – 700X	
MPLFLN10XBDP				
MPLFLN20XBDP				
MPLFLN50XBDP				
MPLAPON50X				
LMPLFLN10XBD				
LMPLFLN20XBD				
LMPLFLN50XBD				



500X	1000X	2000X	5000X	7000X	Working Distance (mm)	NA	Field of View (μm)
					51.7	0.03	19,200 – 2,740
					66.1	0.09	9,100 – 910
140 – 1400X					41.1	0.20	2,740 – 270
					30.0	0.09	9,100 – 910
140 – 1400X					30.0	0.30	2,740 – 270
280 – 2800X					20.0	0.40	1,370 – 140
560 – 5600X					4.5	0.80	690 – 70
					3.5	0.04	17,100 – 2,190
					10.7	0.08	10,200 – 1,100
					12.0	0.15	5,480 – 550
140 – 1400X					6.5	0.25	2,740 – 270
280 – 2800X					3.0	0.40	1,370 – 140
700 – 7000X					1.0	0.75	550 – 55
700 – 7000X					0.35	0.95	550 – 55
140 – 1400X					10.0	0.25	2,740 – 270
280 – 2800X					12.0	0.40	1,370 – 140
700 – 7000X					10.6	0.50	550 – 55

\*The DSX10-SXLOB1, 3, 10X, and DSX10-XLOB3X do not support PO observation.

\*The MPLAPON50X does not support DF and mixed observations.

\*The MPLFLN1.25, 2.5X support BF and OBQ observations.

\*Field of view: At aspect ratio 1:1 diagonal (with factory default value)


### Olympus lens processing system

We created an automatic lens processing system to deliver the highest possible quality optics. As a result, we are now able to process high precision lenses as fine as 1/10,000 mm.



### Olympus advanced engineer development program leads to Yellow Ribbon Medal

In 2018, Olympus was awarded a Yellow Ribbon Medal for developing an advanced method to process high-precision objective lenses up to 2 μm. As part of the program, senior engineers mentored younger engineer in the art and science of lens manufacturing.



# Specifications

## Main Unit Specifications

		DSX10-SZH	DSX10-UZH	
Optical system	Optical system	Telecentric optical system		
	Zoom ratio	10X (motorized)		
	Zoom magnification method	Motorized		
	Calibration	Automatic		
	Lens attachment	Quick-switch, coded lens attachments automatically update magnification and visual field information		
	Maximum total magnification (on monitor)	7,000X		
	Working distance (W.D.)	66.1 - 0.35 mm		
	Accuracy and repeatability (X-Y plane)	Accuracy* <sup>1</sup>	± 3%	
		Repeatability 3 $\sigma_{n-1}$	2%	
Repeatability (Z axis)* <sup>2</sup>	Repeatability $\sigma_{n-1}$	1 $\mu$ m		
Camera	Image sensor	1 / 1.2 inch, 2.35 million pixel color CMOS		
	Cooling	Peltier cooling		
	Frame rate	60 fps (maximum)		
	Normal	1,200 × 1,200 (1:1) / 1,600 × 1,200 (4:3)		
	Fine	Not available	1,200 × 1,200 (1: 1) / 1,600 × 1,200 (4: 3)	
Illumination	Super fine	Not available	3,600 × 3,600 (1: 1) / 4,800 × 3,600 (4: 3)	
	Color light source	LED		
	Lifetime	60,000 h (design value)		
Observation	BF (brightfield)	Standard		
	OBQ (oblique)	Standard		
	DF (darkfield)	Standard LED ring divided into four divisions		
	MIX (brightfield+darkfield)	Standard Simultaneous observation of BF + DF		
	PO (polarization)	Standard		
	DIC (differential interference)	Not available	Standard	
	Contrast up	Standard		
	Depth of focus up function	Not available	Standard	
Focus	Transmitted lighting	Standard* <sup>3</sup>		
	Focusing	Motorized		
	Stroke	101 mm (motorized)		

\*1 Calibration by Olympus or dealer service technician necessary. To guarantee the accuracy of XY, calibration with DSX-CALS-HR (calibration sample) is required.

\*2 When used 20X or higher objective.

\*3 The optional DSX10-ILT is required.

Objective		DSX10-SXLOB	DSX10-XLOB	UIS2
Objective lens	Maximum sample height	50 mm	115 mm	145 mm
	Maximum sample height (free angle observation)	50 mm		
	Parfocal distance	140 mm	75 mm	45 mm
	Lens attachment	Integrated with lens	Available	
	Total magnification	20X - 1,400X	42X - 5,600X	23X* <sup>4</sup> - 7,000X
	Actual F.O.V.	19,200 $\mu$ m - 270 $\mu$ m	9,100 $\mu$ m - 70 $\mu$ m	17,100 $\mu$ m - 50 $\mu$ m
Adaptor	Diffusion adaptor (option)	Available	Not available	
	Eliminate reflection adaptor (option)	Available	Not available	
Lens attachment	Number of objectives that can be attached	Up to 1 piece (attachment is integrated with lens)	Up to 2 pieces	
Objective lens case		Three lens attachments can be stored		

\*4 Total magnification when using MPLFLN1.25X

Stage		DSX10-RMTS	DSX10-MTS	U-SIC4R2
Stage	XY stage: motorized / manual	Motorized (with rotation function)	Motorized	Manual
	XY stroke	Stroke priority mode : 100 mm × 100 mm Rotation priority mode : 50 mm × 50 mm	100 × 100 mm	100 × 105 mm
	Rotation angle	Stroke priority mode : ±20° Rotation priority mode : ±90°	Not available	
	Display rotation angle	GUI	Not available	
	Load-resistance	5 kg (11 lb)		1 kg (2.2 lb)

Frame	DSX-UF	DSX-TF	Display	23 - inch flat panel display
Z-axis stroke	50 mm (manual)		Resolution	1,920 (H) × 1,080 (V)
Tilt observation	Not available	±90°		
Tilt angle display	Not available	GUI		
Tilt angle method	Not available	Manual, fix / release handle		

System Total		Upright frame system	Tilt frame system
Weight (frame, head, motorized stage, display, and console)		43.7 kg (96.3 lb)	46.7 kg (103 lb)
Power consumption		100 - 120V / 220 - 240 V, 1.1 / 0.54A, 50 / 60Hz	

# Customized Solutions

## Expand Your Inspection Capabilities

The DSX1000 digital microscope's precision and ease of use make it a convenient choice for many industrial inspections, and its customization options provide even greater flexibility. Inspections are rarely standard, and a customized DSX1000 microscope can provide the capabilities you need for your application and workflow.

### Beyond Standard

- Larger stages for big and heavy samples
- More space for tall samples without losing image quality
- Added observation modes, such as fluorescence
- And many other customization options



To learn how DSX1000 customized solutions can help you, get in touch:

[www.olympus-ims.com/contact-us](http://www.olympus-ims.com/contact-us)

- OLYMPUS CORPORATION is ISO14001 certified.
- OLYMPUS CORPORATION is ISO9001 certified.
- All company and product names are registered trademarks and/or trademarks of their respective owners.
- Performance characteristics and other values described in this brochure are based on Olympus evaluations as of September 2017, and are subject to change without notice.
- The information including guaranteed accuracy in this brochure is based on the condition set by Olympus. For details, refer to the Instruction Manual.
- Images on the PC monitors are simulated.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

[www.olympus-ims.com](http://www.olympus-ims.com)

**OLYMPUS**<sup>®</sup>

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 2-3-1, Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan