# INSTRUCTIONS

Coded Quintuple Revolving Nosepiece

Coded Quintuple Revolving Nosepiece for BE/DF with Slider Slot for DIC

## U-5RES-ESD / U-D5BDRES-ESD

Connecting the coded revolving nosepiece to the control box for coded units (U-CBS) makes it possible to read out the information in the revolving nosepiece using the U-CBS.


Install the product on the microscope by the following steps.

For the operating procedure, refer to the instruction manuals for the U-CBS (CODED FUNCTION SYSTEM).

Before disposing of this product, be sure to follow the regulations and rules of your local government.

Operating environment

	Temperature	Humidity	Notes
Normal	0 ~ 40°C	30 ~ 90%	
ESD	10 ~ 35°C	40 ~ 70%	It is recommended to use in the clean room.



### 1 Installation of U-CBS


- 1) 4 pieces of rubber feet are included in the coded revolving nosepiece. Stick the rubber feet ① at 4 places.
- 2) Stick the cable holder ② included in the U-CBS.
  - ◎ Be sure that the cable holder is to be stuck at right position as it is illustrated. If it is stuck at connector side for coded revolving nosepiece, cables can not be attached.
- 3) Stick the seal ③ included in the coded revolving nosepiece.
- 4) Place the U-CBS on the table with rubber feet side down.

### 2 Layout of the cable

After installation of the coded revolving nosepiece, stick the cable holders on the microscope and lay out the cable as they are illustrated.

Optical Microscope Accessory


### Mounting the corded Revolving Nosepiece on the STM6

◎ Before mounting the nosepiece, remove the reflected light arm from the STM6 frame.

(As to mount and removal method of the reflected light arm, refer to the instruction manual for the STM6)

- 1) Loosen and remove the screws ② (4 pcs.) of the revolving nosepiece mount on the reflected light arm ① using the Allen screwdriver provided with the microscope frame.
- 2) Turn the revolving nosepiece mount ③ 90 degrees.
- 3) Fix the revolving nosepiece mount by tightening the screws ② (4 pcs. once removed in the above step1) using the Allen screwdriver.
- 4) After mounting the reflected light arm ① on the STM6 frame, fit the coded revolving nosepiece into the mount from left side.

This product is a part of the system product which complies with the CE marking. Please refer to the instruction manual of your system product for the safety instructions related to the CE marking.

In accordance with European Directive on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local EVIDENT distributor in EU for return and/or collection systems available in your country.

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Manufactured by  
**EVIDENT CORPORATION**

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