BIOS/OS settings of PC

BIOS Settings

How to start BIOS:

Turn ON the power of the controller PC, and press the [F10] key while the HP logo is displayed. Depending on the timing to recognize the monitor, the display of the HP logo may not be confirmed. If the OS starts due to the delay in pressing the [F10] key, restart the controller PC and try again.

Setting items			Setting (HP Z2 SFF G9 Workstation) BIOS 01.02.04 / 02.01.02 or later		
S	Security				
	Trusted Execution Technology(TXT)	-	Uncheck		
Advanced					
	Boot Options	Fast Boot	Uncheck ^{*1}		
		Network (PXE) Boot	Uncheck ^{*1}		
	Built-In Device Options	Wake On LAN	Boot to Hard Drive		
	Power Management Options	Extended Idle Power States	Uncheck ^{*1}		
	Remote Management Options	Intel Management Engine (ME)	Uncheck ^{*1}		
		Intel Active Management Technology (AMT)	Uncheck		
		Watchdog Timer	Uncheck (Cannot change)		

*1 Item to change manually (item to be changed from BIOS default)

After confirming that all the above items are set, save the BIOS setting.

How to save BIOS setting:

1. Click [Main] - [Save Changes and Exit], and select [Yes] on the [Save Change?] screen.

2. If the [Secure Boot] screen appears, enter the displayed number and press the [Enter] key.

OS (Windows10) Settings

Before setting the OS, you must activate Windows by connecting the controller PC to the network.

- 1. Connect the controller PC to the network.
- 2. Click [Start] [Settings] [Update & Security] [Activation] [Troubleshoot].
- 3. When [Troubleshoot has completed] appears, click [Activate Windows].
- 4. Make sure that [Windows 10 is Activated with a Digital license] is display, and close all of the setting screen.
- 5. Disconnect the controller PC from the network.

Confirmation items	Setting	
Log on		
When setting up the Windows OS, set the account name and password. After OS setup is completed, make sure that screen lock can be released with [Start] - [Able to logon with the following: Account: "olympus" Password: "olympus"	
BIOS		
Right-click [Start], select [Search], input "SystemInfo", and click [System Information].	BIOS Version in BIOS Settings	
OS		
Edition	Windows10 Pro	
Version	21H2 / 22H2	
Windows setup		
Country or region [Start] - [Settings] - [Time & Language] - [Region]	United States	
Time and date settings [Start] - [Settings] - [Time & Language] - [Region] - [Additional date,time & regional settings] - [change date,time,or number formats]	Match Windows display language (recommended)	
Keyboard layout [Start] - [Settings] - [Time & Language] - [Language] - [Keyboard]	Use language list (recommended)	
TimeZoon Setting [Start] - [Settings] - [Time & Language] - [Data & time] - [Time zone]	Osaka, Sapporo, Tokyo	
[Display] property setting		
Screen saver [Start] - [Settings] - [Personalization] - [Lock screen] - [Screen saver settings]	Screen saver -> None	
	On resume display logon screen -> Off	
TaskBar [Start] - [Settings] - [Personalization] - [Taskbar]	Lock the task bar -> On (No change from default setting)	
	Automatically hide the taskbar in desktop mpde -> On	
	Taskbar location on screen -> Bottom (No change from default setting)	

Confirmation items	Setting
Background Insert the installation medium into the disk drive. Copy the "OLYMPUS_Microscope_2560 x 1600.bmp" in the medium folder (D:\DSX10-BSW\Kitting) into the following folder. C:\Users\Public\Pictures (If your DVD drive is not the D drive, change "D" to the character of your DVD drive.) Right-click the image, select [Set as desktop background].	OLYMPUS OLYMPUS_Microscope_2560 x 1600.bmp
Lock (Sign out) wallpaper Insert the installation medium into the disk drive. Copy "BACKGROUNDDEFAULT.jpg" in the medium folder (D:\ DSX10-BSW\Kitting) into the following folder. C:\Windows\Web\Screen Click [Start] - [Settings] - [Personalization] - [Lock screen]. Click the down arrow button in the background, select [Image] and click [Browse]. Select "BACKGROUNDDEFAULT.jpg" in C:\Windows\ Web\Screen.	BACKGROUNDDEFAULT.jpg
Login account icon Insert the installation medium into the disk drive. Copy "OLYMPUS_icon.bmp" in the medium folder (D:\DSX10-BSW\ Kitting) into the following folder. C:\Users\Public\Pictures Click [Start] - [Settings] - [Accounts]. Click [Browse] under [Create your picture] - [Your picture], select "OLYMPUS_icon.bmp" in C:\Users\Public\Pictures.	OLYMPUS OLYMPUS_icon.bmp
Security setting	
Windows Update [Start] - [Windows System] - [Command Prompt (right-click)] - [More] - [Run as administrator] - Input "gpedit.msc" and press the Enter key - [Local Group Policy Editor] is opened. [Computer Configuration] - [Administrative Templates] - [Windows Components] - [WindowsUpdates]	Configure Automatic Updates: Disable
Microsoft Defender [Start] - [Windows System] - [Command Prompt (right-click)] - [More] - [Run as administrator] - Input "gpedit.msc" and press the Enter key - [Local Group Policy Editor] is opened. [Computer Configuration] - [Administrative Templates] - [Windows Components] - [Microsoft Defender Antivirus] - [Real-time Protection]	Turn off real-time protection: Enabled
Performance setting	
Select a power plan [Start] - [Windows System] - [Control Panel] - [Hardware and Sound] - [Power Options]	High performance
[Change plan settings] - [Turn off the display]	Never
[Change plan settings]] - [Put the computer to sleep]	Never
[Change advanced power settings] - [Hard disk] - [Turn off hard disk after]	"O"
[Change advanced power settings] - [Sleep] - [Sleep after]	Never
[Change advanced power settings] - [Sleep] - [Hibernate after]	Never

Confirmation items		Setting
	[Change advanced power settings] - [Sleep]] - [Allow wake timers]	Disable
	[Change advanced power settings] - [USB settings] - [USB selective suspend setting]	Disable
	[Change advanced power settings] - [PCI Express] - [Link State Power Management]	Off
	[Start] - [Windows System] - [Control Panel] - [Hardware and Sound] - [Power Options] - [Choose what the power buttons do]	Turn on fast startup: Uncheck
	USB hub power settings USB hub power settings [Device Manager] - [Universal Serial Bus controllers] - [USB Root Hub(3.0)] - [Properties] - [Power Management]	Allow the computer to turn off this device to save power: Uncheck
	COM Settings [Device Manager] - [Ports (COM&LPT)] - [Communications Port (COM1)] - [Port Settings] Tab - [Advanced] - [COM Port Number]	COM1
	.NETFramework [Start] (right-click) - [Run] - Input "regedit" and press the Enter key - [Registry Editor] is opened. [Computer] - [HKEY_LOCAL_MACHINE] - [SOFTWARE] - [Microsoft] - [NET Framework Setup] - [NDP] - [v4] - [Full] -> Check value of REG_DWORD described in [Release].	.NET4.8 (value of Release DWORD: 528449)
G	aphic settings	
	[Start] - [NVIDIA Control Panel] - [Adjust Image Settings with Preview]	[Use my preference emphasizing]: Performance
	[Start] - [NVIDIA Control Panel] - [View system topology]	Ver.473.47
Ad	dditional user (LOG IN user)	
	Computer Name [Start] - [Windows System] - [Control Panel] - [System and Security] - [System] - [Change settings] - [Computer Name] - [Change] Changing computer name, select [Restart later].	Computer Name: OLYMPUS-PC
	Interactive logon [Start] (right-click) - ([Search]) - enter "local" - click [Local Security Policy] - [Local Policy] - double click [Security Options]	Interactive logon: Don't display last signed-in: Disabled (No change from default)
	Auto login Insert the installation medium into the disk drive. Double-click AutoLogin.reg in the following folder in the installation medium. D:\DSX10-BSW\Kitting Set according to the displayed registry editor screen.	Automatic logon to the operating system
In	stallation software	
	Pre-installed app	PDF browser (Microsoft Edge)

OS (Windows11) Settings

Before setting the OS, you must activate Windows by connecting the controller PC to the network.

- 1. Connect the controller PC to the network.
- 2. Click [Start] [Settings] [Update & Security] [Activation] [Troubleshoot].
- 3. When [Troubleshoot has completed] appears, click [Activate Windows].
- 4. Make sure that [Windows 11 is Activated with a Digital license] is display, and close all of the setting screen.
- 5. Disconnect the controller PC from the network.

Confirmation items	Setting	
BIOS		
[Start] (right-click) - [Search] - Input "SystemInfo" Click [System Information].	BIOS Version or later in BIOS Settings	
OS		
edition [Start] - [Setting] - [System] - [About]	Windows 11 Pro	
version [Start] - [Setting] - [System] - [About]	21H2 / 22H2	
Windows setup		
Region [Start] - [Settings] - [Time & language] - [language & reginon] - [Windows display language]	United States	
Time & language [Start] - [Settings] - [Time & language] - [language & reginon] - [Regional format] or [Start] - [Time & language] - [language & reginon] - [Administrative languages settings] - [Formats] tab	Match Windows display language (recommended)	
Keyboard layout [Start] - [Settings] - [Time & language] - [Typing] - [Advanced keyboard settings] - [Override for default input method]	Use language list (recommended)	
Setting time zone [Start] - [Settings] - [Time & language] - [Data and Time] - [Time zone]	Local time	
[Display] property setting		
Screen saver	Screen saver: None	
[Start] - [Settings] - [Personalization] - [Lock screen] - [screen saver]	On resume, display logon screen: OFF	
Taskbar	Automatically hide the taskbar: ON	
[Start] - [Settings] - [Personalization] - [Taskbar] - [Taskbar behaviors]	Taskbar alignment: Center (No change from default)	

С	onfi	rmation items	Setting
	Ba Ins Co me C: (If y yo Rig	ackground sert the installation medium into the disk drive. opy the "OLYMPUS_Microscope_2560 x 1600.bmp" in the edium folder (D:\DSX10-BSW\Kitting) into the following folder. \Users\Public\Pictures your DVD drive is not the D drive, change "D" to the character of ur DVD drive.) ght-click the image, select [Set as desktop background].	OLYMPUS OLYMPUS_Microscope_2560x1600.bmp
	LC Ins Cc Ci Cli Cli an Se Sc	ock (Sign out) wallpaper sert the installation medium into the disk drive. opy "BACKGROUNDDEFAULT.jpg" in the medium folder (D:\ SX10-BSW\Kitting) into the following folder. \Windows\Web\Screen ick [Start] - [Settings] - [Personalization] - [Lock screen]. ick the down arrow button in the background, select [Image] id click [Browse]. elect "BACKGROUNDDEFAULT.jpg" in C:\Windows\Web\ creen.	BACKGROUNDDEFAULT.jpg
	Lo Ins Co BS C: Cli Cli Cli	gin account icon sert the installation medium into the disk drive. by "OLYMPUS_icon.bmp" in the medium folder (D:\DSX10- SW\Kitting) into the following folder. \Users\Public\Pictures ick [Start] - [Settings] - [Accounts]. ick [Browse] under [Create your picture] - [Your picture], select 'LYMPUS_icon.bmp" in C:\Users\Public\Pictures.	OLYMPUS OLYMPUS_icon.bmp
Security setting		rity setting	
	Se [St ad - [[[Co Co ex	etting Windows Update [art] - [All apps] - [Windows Tools] - [Command Prompt (Run as Iministrator)] - Input "gpedit.msc" and press the Enter key Local Computer Policyr] is opened. omputer Configuration] - [Administrative Templates] - [Windows omponents] - [Windows Update] - [Manage end user perience]	Configure Automatic Updates: Disable
	Wi [St ad - [L [Co Co Pro	indows Defender art] - [All apps] - [Windows Tools] - [Command Prompt (Run as Iministrator)] - Input "gpedit.msc" and press the Enter key Local Computer Policyr] is opened. omputer Configuration] - [Administrative Templates] - [Windows omponents] - [Microsoft Defender Antivirus] - [Teal-time otection]	Turn off real-time protection:Enabled
Performance setting			
	Se [St an	elect power options art] - [All apps] - [Windows Tools] - [Control Panel] - [Hardware ad sound] - [Power Options] - ([Show additional plans])	High performance
		[Change plan settings] - [Turn off the display]	Never
		[Change plan settings] - [put the computer to sleep]	Never
		[Change advanced power settings] - [Hard disk] - [Turn off hard disk after]	"O"

Confirmation items		rmation items	Setting
		[Change advanced power settings] - [Sleep] - [Hibernate after]	Never
		[Change advanced power settings] - [Sleep] - [Allow wake timers]	Disable
		[Change advanced power settings] - [USB settings] - [USB selective suspend setting]	Disable
		[Change advanced power settings] - [PCI Express] - [Link State Power Management]	Off
	Sh [Si ar do	nutdown settings tart] - [All apps] - [Windows Tools] - [Control Panel] - [Hardware nd sound] - [Power Options] - [Choose what the power buttons o] - [Change settings that are currently unavailable]	[Turn on fast startup] checkbox: no check
	[Si Co M	tart (right-click)] - [Device Manager] - [Universal Serial Bus ontrollers] - [USB Root Hub(USB3.0)] - [Properties] - [Power anagement]	Allow the computer to turn off this device to save power: no check
	[Si [C [C	tart (right-click)] - [Device Manager] - [Ports (COM & LPT)] - communications Port (COM1)] - [Port Settings] - [Advanced] - cOM Port Number]	COM1
	.N [Si - [[C - [->	ETFramework tarl] (right-click) - [Run] - Input "regedit" and press the Enter key Registry Editorr] is opened. omputer] - [HKEY_LOCAL_MACHINE] - [SOFTWARE] - [Microsoft] NET Framework Setup] - [NDP] - [v4] - [Full] Check value of [REG_DWORD] described in [Release].	HP Z2 SFF G9 Workstation: .NET4.8 (value of Release: 528449 / 533320)
Graphic settings			
	[Start] - [All apps] - [NVIDIA Control Panel] - [Adjust image settings with preview]		Use my preference emphasizing: Performance
	[S	tart] - [All apps] - [NVIDIA Control Panel] - [View system topology]	Driver version: Ver.473.47
Additional user (LOG IN user)			
	Computer name: [Start] - [Setting] - [System] - [About] - [Rename this PC] Changing computer name, select [Restart later].		Computer name: Serial Number
	Int [St Pc	teractive logon: tart] (right-click) - ([Search])- enter "local" - click [Local Security plicy] - [Local Policy] - double click [Security Options]	Interactive logon: Don't display last signed-in: Disabled (No change from default)
Installation software			
	Pr [Si	e-installed application tart]	The following is pre-installed. PDF browser (Microsoft Edge)

EVIDENT CORPORATION