BXC-CBB Main control box for BXC-CBB system.



SPECIFICATIONS

STECHTCATIONS				
Global item code / Local code for USA		N6387900 / BXC-CBB		
Predecessor product information		None		
Product Category		Motorized Units		
Rating	Control box body	Input: DC 24 V 2.5 A(Max)		
	AC adapter	Input: AC 100-240 V 50-60 Hz 1.4 A(Max) Output: DC 24 V 2.71 A(Max)		
Cooling system		Natural cooling		
Cooling system Fixing method		 Installation on the floor with rubber legs. Can be fixed on the floor, wall, ceiling, etc. No space necessary around the unit. The recommended screw length is 8 mm or more. Recommended to remove rubber legs when fixing with screws. 		

BXC-CBB

RS232C communication control I/F	[Communication form] RS-232C [Connector] • D-Sub9pin DCE array (front panel) • Clamping screw for fitting: #4-40 UNC • Name of connec • tor on back panel: "RS-232C" [Required cables] • Straight connection type with D-Sub9pin(female) - D-Sub9pin(female) [Settings for communication]* Fixed		in(female)
	Baudrate	19200 [bps]	
	Data bit	8 [bits]	
	Parity	even	
	Stop bit	1 or 2 [bits] (No need to switch settings)*	
	Terminator	CR+LF	
	Flow control	None	
	* Setting of Sto StopBit is fixed t	pBit when communicating from Host-PC to co "2" when communicating from BXC-CBB t	BXC-CBB. o Host-PC.
Weight	About 0.4 kg (excluding AC adapter) *AC adapter: About 0.3 kg		
Country of Origin	Japan		

SYSTEM ENVIRONMENT LIMITATIONS FOR OPERATION

Environment	Indoor use
Temperature T0 [°C]	5 – 40
Humidity RH0 [%]	30 - 85

RESTRICTIONS


The total length of cables used (BXC-LCBL1M, BXC-LCBL3M, and BXC-LCBL6M) is up to 13 m.


PACKING LIST


Item	Quantity
Control box	1
AC adapter	1

MAIN DIMENSION

AC adapter: 53(W) x 38(H) x 115(D) mm (excluding protrusions)








(Unit: mm)

BXC-CBB

SYSTEM DIAGRAM



*1 The BXC-FSU can be connected directly to the BXC-CBB to use auto focus with a manual illuminator. In this case, please note that the BXC-CBB cannot control the BX3M-LEDR without the BXC-RLI. If you want to use the BX3M-LEDR with a manual illuminator, please also use the BXC-CBRML or use the BX3M-PSLED, which is a local control.

*2 When using BXC-FSU, DIC observation and simple polarizing observation are not possible even if it is manually operated.

· Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.